

DATA SHEET SOLO Plus – Personal Bleed Control Pack

SKU Number **SS-6010150040**
 EAN13 Barcode **5060483286402**

Date of Publication **09-01-2025**
 Commodity Code **30065000.00**



| | |
|----------------------------|--|
| <p>Product Description</p> | <p>SOLO Plus – Personal Bleed Control Pack</p> <p>Effective haemorrhage control requires the right equipment readily available to ensure first responders can act swiftly.</p> <p>The STEIN SOLO Plus Personal Bleed Control Pack is designed to deliver essential, high-quality medical equipment for managing a life-threatening haemorrhage. This kit enables you to provide potentially lifesaving care whilst awaiting Emergency Services.</p> <p>Contained in a robust, weather-resistant pouch, the kit features a sealed zipper for added protection. It is equipped with a Velcro belt loop, a securing press stud, and a D-ring, allowing it to attach securely to any climbing harness for convenient access.</p> |
|----------------------------|--|

| | |
|-----------------|--|
| Contents | <p>1 x 4" Emergency FCP Bandage - Military</p> <p>1 x WoundClot Haemostatic Gauze Dressing 10 x 10</p> <p>1 x Eye/Wound Solution (20ml Ampoule)</p> <p>1 x Disposable Medical Gloves</p> <p><u>Plus Version Includes</u></p> <p>1 x SWAT-T Rescue Orange Tourniquet</p> |
|-----------------|--|

| | |
|-------------------------|---|
| Safety Standards | <ul style="list-style-type: none"> All medical products used in the STEIN Kits are CE Marked and conform with the relevant European Union harmonisation legislation. |
|-------------------------|---|

| | |
|-----------------------|--|
| Product Weight | <ul style="list-style-type: none"> 284g |
| Size | <ul style="list-style-type: none"> 100cm x 140cm x 50cm (approx.) |

| | |
|-----------------|--|
| Storage | <ul style="list-style-type: none"> Keep stored and transported at ambient temperature (0°C to + 30°C) Do not store in direct/high heat or sunlight as this may affect the internal components of the kit |
| Lifetime | <ul style="list-style-type: none"> Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable. |
| Cleaning | <ul style="list-style-type: none"> Clean the storage bag using warm soapy water and soft cloth |
| Disposal | <ul style="list-style-type: none"> Recycling is preferred where possible. Look out for the recycling symbol for material categories |

| | |
|-------------------------------|--|
| Replacement Components | <ul style="list-style-type: none">• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents• Refer to the STEIN Website for replacement items www.steinworldwide.com |
| Additional Information | <ul style="list-style-type: none">• In extreme cases of haemorrhage further products may be required over and above the items in these kits.• This kit does not contain products for the treatment of minor injuries• This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019 |
| Training | <ul style="list-style-type: none">• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use.• Training Providers can be recommended upon request |
| Trademarks | <ul style="list-style-type: none">• STEIN® is a registered trademark |