

DATA SHEET SOVOS Platinum FORCE

SKU Number **SV-95KBC00050**
 EAN13 Barcode **5060660858125**

Date of Publication **09-01-2025**
 Commodity Code **30065000.00**



| | |
|-----------------------------------|---|
| <p>Product Description</p> | <p>SOVOS Platinum FORCE Bleed Control Trauma Kit</p> <p>When controlling severe haemorrhages, having a kit equipped with the right medical supplies is essential. User-friendly Bleed Control Trauma Kits are vital to the effectiveness of first responders in critical situations.</p> <p>The SOVOS Platinum FORCE Bleed Control Trauma Kit is designed to handle life-threatening haemorrhages, lacerations, Burns, and minor injuries providing high-quality medical supplies to meet urgent needs. Packed in a compact, lightweight bag, the kit can be secured to a harness for immediate and rapid access. The kit is supplied with a load rated tether and connector.</p> <p>Packed in a compact yet spacious bag, this kit offers increased capacity to accommodate additional critical materials essential for first responders. The fast-access, single-pull, full-length front opening flap ensures immediate visibility and access to all contents.</p> <p>A removable, double-sided central storage panel enhances organization, featuring mesh-fronted zip pockets on one side and elasticated loops on the reverse.</p> <p>Equipped with anti-slip leg straps, a removable waistband, and PAL webbing, the pack can be securely attached to the body or integrated with other gear, ensuring ease of use and reliable positioning in high-pressure situations.</p> <p>The term "Platinum 10" highlights the critical first 10 minutes after a traumatic injury or medical emergency. This window of time is vital for swift and efficient action by first responders, emergency personnel and potentially the casualty themselves to improve survival rates and minimise long-term harm through controlling major bleeds and stabilising injuries.</p> <p>The SOVOS Platinum kits enable trained first aiders, advanced medical professionals and first responders to deliver lifesaving care during emergencies while awaiting the arrival of the emergency service team.</p> |
|-----------------------------------|---|

| | | | |
|-----------------|---|----------------------|--------------------------------------|
| Contents | 1 | x | Platinum FORCE Storage Bag |
| | 1 | x | Load Rated Tether with Connector |
| | 1 | x | SAM XT Tourniquet |
| | 1 | x | 6" FCP Military Emergency Bandage |
| | 1 | x | ChitoSAM Haemostatic Gauze 6' Z-Fold |
| | 1 | x | Medical Gauze No14 - Wound Dressing |
| | 1 | x | Large First Aid Dressing 18cm x 18cm |
| | 1 | x | Burns Dressing 10cm x 10cm |
| | 1 | x | Emergency Foil Thermal Blanket |
| | 1 | x | Sterile Saline Spray (100ml) |
| | 1 | x | 20ml Ampoule Sterile Wound Wash |
| | 1 | x | Resuscitation Face Shield |
| | 1 | x | Pack of 10 Plasters |
| | 1 | x | Pack Sterile Skin Closure Strips (5) |
| | 1 | x | High Performance Medical Shears |
| | 1 | x | Medical Marker Pen |
| | 2 | x | Disposable Medical Gloves |
| | 2 | x | De-Contamination Wipes 30cm x 20cm |
| | 2 | x | Sterile Wound Cleansing Wipes |
| | 1 | x | Emergency Whistle |
| 1 | x | Emergency Glow Stick | |

| | |
|-------------------------|--|
| Safety Standards | <ul style="list-style-type: none"> All medical products used in the SOVOS Platinum Kits are CE Marked and conform with the relevant European Union harmonisation legislation. |
|-------------------------|--|

| | |
|-----------------------|---|
| Product Weight | <ul style="list-style-type: none"> 1218g 1350g (With leg straps fitted) |
| Size | <ul style="list-style-type: none"> 10cm x 21cm x 20cm (approx.) |

| | |
|-----------------|--|
| Storage | <ul style="list-style-type: none"> Keep stored and transported at ambient temperature (0°C to + 30°C) Do not store in direct/high heat or sunlight as this may affect the internal components of the kit |
| Lifetime | <ul style="list-style-type: none"> Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable. |
| Cleaning | <ul style="list-style-type: none"> Clean the storage bag using warm soapy water and soft cloth |
| Disposal | <ul style="list-style-type: none"> Recycling is preferred where possible. Look out for the recycling symbol for material categories |

| | |
|-------------------------------|---|
| Replacement Components | <ul style="list-style-type: none">• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents• Refer to the SOVOS Website for replacement items www.sovoshelmets.com |
| Additional Information | <ul style="list-style-type: none">• In extreme cases of haemorrhage further products may be required over and above the items in these kits.• This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019 |
| Training | <ul style="list-style-type: none">• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use.• Training Providers can be recommended upon request |
| Trademarks | <ul style="list-style-type: none">• SOVOS® is a registered trademark |